


## General Description of Vitreous Carbon Discs and Graphite Discs



The standard SEM stubs used in an SEM are nearly all made from aluminum and sometimes brass. This is fine for many applications if the sample is large enough and can be mounted directly on the stub. However, for small samples, powders, particles in solutions and fibers the finish and the material of the stub can interfere with imaging or X-ray analysis.

### Carbon (graphite) or vitreous carbon discs to prevent interference.

To prevent interference of the stub surface you can place a carbon (graphite) disc or a vitreous carbon disc between the stub and the specimen and mount the specimen directly on these discs. Both discs varieties (also called planchets) form a protective shield.

### Differences between vitreous carbon discs and carbon (graphite) discs.

Vitreous carbon is superior in hardness, surface finish, mechanical and chemical stability but it is more expensive than carbon discs. The table below will give guidance to which material is best for the application at hand. Generally speaking, vitreous carbon is the better choice, but for less demanding applications carbon (graphite) discs will do just fine.

	Vitreous Carbon Disc	Carbon (Graphite) Disc
Smooth surface	V	X
Hard Surface	V	X
Easy to clean and re-use	V	O
SE Imaging	V	X
BSE Imaging	V	V
X-ray Microanalysis	V	V
Spectral Analysis	V	O
Aqueous Solutions	V	X
Pricing	O	V
High Purity	V	V
Conductivity	V	V
Packaging	Each	Pkg/10

## Vitreous Carbon Discs

Vitreous carbon, glassy carbon or glass-like carbon has many application in research, science, metallurgy and chemical analysis. For microscopy related techniques this material has a number of desirable properties which makes it an ideal substrate material:

- Vitreous carbon combines the glassy, ceramic properties of carbon with those of graphite
- It is a conductive material
- It is a high purity, hard, non-porous, corrosion resistant material
- Impermeable to gasses and / or liquids
- Good thermal stability and great resistance against thermal shock
- Smooth surface
- Less hydrophobic than graphite

For microscopy and analytical applications it is important to know that the smooth surface gives rise to low electron signals, both for secondary electrons (SE) and backscattered electrons (BSE). The high purity carbon material contributes only C to the X-ray spectrum which makes vitreous carbon ideal for EDX and WDX investigations of small samples, powders, particulates, and fibers. Vitreous carbon acts less hydrophobic than graphite and is more suitable to prepare particles in aqueous solution where the solution will spread more evenly over the surface.

Physical Properties of HG Grade Vitreous Carbon	
Maximum Temperature (vacuum or inert gas), <sup>o</sup> C	3000
Density, g/cm <sup>3</sup>	1.42
Hardness, Vickers, HV	230
Permeability Coefficient, cm <sup>2</sup> /s	10 <sup>-9</sup>
Open porosity, %	0
Young's Modulus, GPa	35
Compressive Strength, Gpa	0.48
Flexural Strength, Gpa	0.26
Thermal Conductivity (30 <sup>o</sup> C) W/km	6.3
Electrical Resistance (30 <sup>o</sup> C), Ωμm	45
Thermal Expansion Coefficient (20-200 <sup>o</sup> C), /K	2.6x10 <sup>-6</sup>

### Surface finish of Vitreous Carbon Discs

The Vitreous carbon discs are lapped on one side to provide a smooth surface, typically 0.05 – 0.01μm. They can be used many times after cleaning and lapping. For lapping we suggest to use either diamond lapping film with water or SiC in water suspension.

### **Purity, handling and sizes of the vitreous carbon discs.**

The high purity vitreous carbon material has less than 30ppm impurities in total. Typical impurities for vitreous carbon can be:

- Ca 10ppm
- Si 14ppm
- Al, Fe, K, Na, Ni and Sn all less than 1ppm
- Ba, Be, Bi, Cd, Co , Cr, Cu, Mg, Mo, Sr, Ti, V, W, Zn and Zr all less than 0.1ppm

The sizes available are 10, 12.7, 19.0 25.4 and 32mm diameters which fit on a the most popular pin stubs and cylinder mounts. The discs have a thickness of either 2 or 3mm and can be easily handled. For SEM application we suggest to use conductive strong carbon glue (preferable) or silver glue to mount the vitreous carbon discs on the SEM stubs.

### **Graphite Discs**

The graphite planchets or discs consist of high purity, graphitized soft spectral grade carbon. This material is softer and is less strong than the vitreous carbon materials. When handled with care is still an excellent substrate for SEM imaging and X-ray micro-analysis. The carbon discs exhibit a low contribution for BSE imaging and a pure background with only the contribution of carbon for EDX and/or WDX analysis. The standard finish is not a smooth as with the vitreous carbon discs and is visible in SE imaging. Carbon discs are also more porous and more hydrophobic than vitreous carbon discs.

<b>Physical Properties of Spectral Grade High Purity Carbon</b>	
Density, g/cm <sup>3</sup>	1.6
Hardness, Shore	33
Porosity, %	29
Young's Modulus, Mpa	22
Compressive Strength, Mpa	48
Flexural Strength, Mpa	27
Electrical Resistance, Ωμm	65
Thermal Expansion, /K	1.9x10 <sup>-6</sup>

### **Surface finish of graphite discs**

The standard ground finish on one side is 0.80μm. The finish can be easily improved by lapping the discs on fine 10um grinding or lapping paper with water as lubricant; the surface roughness of the soft graphite material can be improved to 0.05μm or better.

### **Purity, handling and sizes of graphite discs**

The high purity carbon used for the carbon planchets or discs only has impurities of less than 2ppm with 1ppm or less for each single element. The impurities can be B, Mg, Al, Si, Ca or Fe.

The sizes available are 10, 12.7, 25.4 and 32mm diameters, all with a thickness of 1.6mm. These 4 sizes can be used on all popular SEM stubs and mounts. Since the material is relatively soft, we strongly suggest to glue the planchets on SEM stubs or mounts with a conductive carbon or silver glue. We would like to discourage the use of SEM stubs fully made of carbon; the material is soft so the pin can easily break off and during handling the soft carbon can easily shed material which can potentially cause contamination. Carbon discs glued on SEM stubs combine the advantages of the low contributing carbon background, low costs, easy handling and storage with the sturdy aluminum base. When the top of the specimen stub is completely covered by the graphite discs there is no imaging or analytical contribution of the supporting SEM stub