

## General Description of Preparation Pods and SEM Stub Stands

The Preparation pods and SEM stub stands are designed for secure handling of SEM stubs and to avoid confusion when preparing multiple SEM stubs. For the pin stubs a preparation pod is a necessity to hold the pin stub in the correct position. For the JEOL and Hitachi cylinder stubs the SEM stub stands provide a way to handle and transport multiple stubs. The SEM stub stands are also excellent for storing SEM stubs and or SEM holders in a desiccator or other storage location.

There are two type of preparation pods available. They are made from hard anodized aluminum:

- Preparation pods horizontal stand for 11 standard pin stubs



- Multi-angle preparation pods stand for horizontal, 45° and 90° pin stubs



The stub stands are available for all types of SEM stubs. Made from clean aluminum:



- Multi stub preparation stand for 23 pin stubs
- Multi stub preparation stand for 38 pin stubs
- Multi stub preparation stand for 28 JEOL 9.5mm stubs
- Multi stub preparation stand for 28 JEOL 12.2mm stubs

- Multi stub preparation stand for 28 Hitachi  $\varnothing$ 15mm stubs
- Multi stub preparation stand for 8 JEOL or Hitachi  $\varnothing$ 25mm stubs
- Multi stub preparation stand for 6 JEOL or Hitachi  $\varnothing$ 32mm stubs